


# GeneProof®



## GeneProof COVID-19 Solution



Thanks to a constant process of innovation, we are gradually creating **GeneProof COVID-19 Solution** - this is a series of PCR tests, which aims to offer a comprehensive solution in the field of detection of SARS-CoV-2.



PCR tests for **different types of use** - diagnostics, screening, genotyping.



Unique specificity and sensitivity - **excellent QCMD panel results.**



W.H.O. based design - **avoids fast negative results** caused by virus mutations.



Faster examination with **FAST PCR** profile - **less than 1 hour.**



**Easy-to-Use concept** - single tube Ready-to-Use Master Mix.



Wide range of verified medium types and specimens - **swab, saliva, transport medium.**



**Direct detection** possibility.



**View our COVID-19 Solution and individual products on our website or get a quick overview in our COVID-19 portfolio**



[www.geneproof.com](http://www.geneproof.com)

ISO 13485

GUARANTEED CONTROL OVER THE COMPLETE PROCESS OF DEVELOPMENT, MANUFACTURING AND DISTRIBUTION OF ALL THE OFFERED GENEPROOF PRODUCTS

# GeneProof COVID-19 Solution

GeneProof PCR kits	Diagnostics	Screening	Typing	Direct detection	
				Bi-CoV™ set	Other media
<b>GeneProof SARS-CoV-2 PCR Kit</b> diagnostic test with confirmation in an independent channel	<i>RdRp/E</i> gene, <i>N</i> gene, Internal Control, 3 channel			✓	
<b>GeneProof SARS-CoV-2 Screening PCR Kit</b> screening test with triple-target detection		<i>RdRp/E/N</i> gene, Internal Control, 2 channel		✓	✓
<b>GeneProof SARS-CoV-2 RdRp Confirmatory PCR Kit</b> affordable solution for very fast and simple screening		<i>RdRp</i> gene, Internal Control, 2 channel		✓	✓
<b>GeneProof SARS-CoV-2 Vaccine Escape PCR Kit</b> fast simultaneous SARS-CoV-2 detection and typing of the most dangerous vaccine escape mutations		<i>RdRp</i> gene, <i>S</i> gene, Internal Control, 4 channel	E484K, N501Y mutations		
<b>GeneProof SARS-CoV-2 South African/UK Diagnostic PCR Kit</b> fast simultaneous SARS-CoV-2 detection and typing of variants		<i>RdRp</i> gene, <i>S</i> gene, Internal Control, 4 channel	South African/UK	✓	
<b>GeneProof SARS-CoV-2 South African/UK Typing PCR Kit</b> fast genotyping of SARS-CoV-2 variants			South African/UK, 3 channel		
<b>GeneProof SARS-CoV-2 South African/UK/BR Typing PCR Kit</b> fast genotyping of SARS-CoV-2 variants			South African/UK/BR, 4 channel		

## NON-INVASIVE AND COMFORTABLE SAMPLING

**Direct detection Bi-CoV™ set** – screening solution for massive testing of public based on the non-invasive and comfortable sampling from the front nose or mouth (saliva sample) using the sampling kit Bi-CoV™. You do not need to isolate the viral RNA, therefore the procedure is significantly shorter, easier, and cheaper.



COMPATIBLE WITH A WIDE  
RANGE OF REAL-TIME PCR  
DEVICES



CERTIFIED  
DIAGNOSTIC TESTS